






The ExtremiFix cannulated screw system offers 2.0, 2.4, 3.0 and 4.0mm cannulated screws in both headed and headless screw options to treat a variety of indications in the upper extremity. The ExtremiFix system is carefully engineered to encompass the following features and benefits:

- Titanium screws designed to provide superior bone purchase and excellent compression
- Self drilling and self tapping screws to facilitate the surgical technique
- Reverse cutting flutes to ease screw removal
- · Low profile screw heads reduce soft tissue irritation
- Headless screws to facilitate countersinking, and feature a small trailing thread diameter
- Variety of sizes, from 2.0 mm for hand fractures and fusions to 4.0 mm for elbow injuries
- 1 system with 2 types of screws, allows for intra-operative decisions
- Tapered driver stems to facilitate percutaneous implantation
- 1mm screw implant size options

# **Implant Block Options**








3.0/4.0mm BLOCK



2.0-3.0mm Cannulated BLOCK



## Implant Size Options

#### **CANNULATED HEADED IMPLANTS**

319-20xx 2.0mm x 6mm-14mm (1mm increments) 2.0mm x 16mm-42mm (2mm increments) 2.4mm x 6mm-14mm (1mm increments) 2.4mm x 16mm-50mm (2mm increments) 319-30xx 3.0mm x 10mm-40mm (2mm increments)

319-40xx 4.0mm x 12mm-52mm (2mm increments)

#### **CANNULATED HEADLESS IMPLANTS**

317-20xx 2.0mm x 10mm-14mm (1mm increments) 2.0mm x 16mm-42mm (2mm increments) 317-24xx 2.4mm x 10mm-14mm (1mm increments) 2.4mm x 16mm-50mm (2mm increments) 317-30xx 3.0mm x 10mm-14mm (1mm increments) 3.0mm x 16mm-40mm (2mm increments) 317-40xx 4.0mm x 12mm-52mm (2mm increments)

### OsteoMed Products



ExtremiLOCK Ankle Plating System



ExtremiLOCK Ankle Fusion System



ExtremiLOCK Foot Plating System



Large Cannulated Screws



**HPS Hand Plating System** 



Hand Fusion



Pet Device



CMC Fusion Plate



Inion



OsteoVation EX



OsteoVation QWIK



OsteoVation B-TCP



Allograft Offering