


TRAUSON®

Universal Retrograde Nail System

Brochure



Products you can trust at a price you can afford

Universal Retrograde Nail System

Indications

- The Universal Retrograde Nail System is intended for the fixation of fractures where flexibility of the implant is desired. This includes fractures of: Lower extremity diaphyseal fractures of children and smallstatured adults, upper extremity diaphyseal fractures in both, adults and children as well as some metaphyseal fractures, such as radial neck, proximal humerus and supracondylar humerus fractures. Please refer to the surgical technique for a complete list of indications and contraindications.

Features & Benefits

- Minimal invasive technique and biomechanical stability allow early functional exercise
- Type II anodized Titanium Alloy for enhanced biomechanical and biomedical performance
- Curved nail tip allows maximal guidance without blocking in the intramedullary canal
- Tapered tip and concave inner nail tip for anchoring in metaphyseal bone and capturing fracture fragments

Implants

Universal Retrograde Nail

P/N	Dia	Length	Material
46116204E	Ø2	440	TC4
46116254E	Ø2.5	440	TC4
46116304E	Ø3	440	TC4
46116354E	Ø3.5	440	TC4
46116404E	Ø4	440	TC4

End Cap

P/N
46119001E

Instruments

P/N	Description	Quantity
12600100E	Straight Impactor	1
12600200E	Hammer	1
12600300E	Wrench, Ø4.5	1
12600400E	Pointed Awl	1
12600500E	Pin Holding Clamp	1
12600600E	Insertor	1
12600700E	Three-cornered Wrench	1
12600800E	Cutter	1
12600900E	Oblique Impactor	1
12601000E	Drill Bit, Ø2.7	1
12601200E	Spanning Wrench	1
12601300E	Instrument Case	1

TRAUSON®

Trauson(China) Medical Instrument Co., Ltd

Add: No.9 Long Men Road, Wujin Hi-Tech Industrail Zone, 213164, Changzhou City, Jiangsu Province. P.R.C.
Tel: 86-519-86387075 Fax: 86-519-86387279 Email: info@trauson.com http://www.trauson.com