

3.0 Instruments Ordering Information

P/N	Description	Quantity
14900100E	Ø1.1 Threaded Guide Wire	5
14900200E	Ø1.1 Guide Wire	5
14901201E	Ø2.0 Drill Bit for Quick-connect Handle, Cannulated	2
14900301E	Ø2.0 Drill Bit, Cannulated	2
14901400E	Ø2.0 Drill Bit for Quick-connect Handle	1
14901500E	Ø2.0 Drill Bit	1
14301700E	Handle with Quick Coupling	1
14900500E	Double Drill Guide Ø2.0/Ø1.1	1
14900401E	Stardrive Screwdriver Shaft, Cannulated	1
14900601E	Stardrive Screwdriver Shaft	1
14900700E	Screw Holding Sleeve	1
14900800E	Ø1.1 Cleaning Wire	1
14900900E	Depth Gauge	1
14901301E	Countersink, Cannulated	1
10901600E	Screw Holder	1
10700200E	Periosteal Elevator, Small	1
10700300E	Retractor, Small	1
14901000E	Screw Box	1
14901100E	Outer Case	1

4.5 Instruments Ordering Information

P/N	Description	Quantity
10902001E	Drill Bit, Ø1.5	2
10903100E	Threaded Guide Wire, Ø1.5	3
10902200E	Drill Guide, Ø3.2/Ø1.6	1
10901133E	Cannulated Drill bit	2
10902300E	Cannulated Tap	2
10902400E	Drill Guide, Ø7.0/Ø3.2	1
10902500E	Protective Sleeve	1
10902600E	Depth Gauge	1
10902700E	Cannulated Spot-facer	1
10902800E	Cannulated screwdriver with holding Sleeve	1
10902900E	Parallel Drill Guide	1
10900400E	Canal Cleaner	1
15801400E	Hexagonal Screwdriver, SW3.5	1
10904000E	Extractor for Hexagonal Screw, SW3.5	1
10904400E	Screw Box	1
10904100E	Outer Case	1

4.0 Instruments Ordering Information

P/N	Description	Quantity
15000100E	Threaded Guide Wire Ø1.25	5
15000200E	Guide Wire Ø1.25	5
15000301E	Ø2.7 Drill Bit for Quick-connect Handle, Cannulated	2
15000401E	Ø2.7 Drill Bit, Cannulated	2
15000501E	Ø5.6 Countersink, Cannulated	1
15000601E	Ø4.0 Tap, Cannulated	1
15000900E	Double Drill Guide Ø2.7/Ø1.25	1
15000701E	Ø2.8 Stardrive Screwdriver, Cannulated	1
15001500E	Screw Holding Sleeve	1
15001000E	Depth Gauge	1
15000800E	Ø1.25 Cleaning Wire	1
10901600E	Screw Holder	1
14301700E	Handle with Quick Coupling	1
15001401E	Ø 2.8 Stardrive Screwdriver	1
15001600E	Ø2.7 Drill Bit	1
15001300E	Ø2.7 Drill Bit for Quick-connect Handle	1
15001100E	Screw Box	1
15001200E	Outer Case	1

7.3 Instruments Ordering Information

P/N	Description	Quantity
10900700E	Parallel Drill Guide- I	1
10900800E	Parallel Drill Guide- II	1
10900500E	Protection Sleeve, Ø11/Ø8	1
10900600E	Guide Wire Sleeve, Ø8/Ø2.6	1
10903200E	Threaded Guide Wire, Ø2.3	3
10704500E	Drill Bit, Ø2.5 mm	1
10900900E	Depth Gauge	1
10901006E	Drill Bit, Cannulated, Ø5	1
10901200E	Tap, HB7.3	1
10901300E	Cannulated Spot-facer, HB7.3	1
10901500E	SW4.0 Hex Screwdriver, Cannulated	1
10901600E	Screw Holder	1
10900300E	Canal Cleaner, Ø2×260	1
10902100E	Screwdriver, SW4	1
10903900E	Extractor for Hex Screw, SW4	1
10904600E	Screw Box	1
10904200E	Outer Case	1

Cannulated Screw

Brochure




Cannulated Screw

Indications

- Cannulated screws are indicated for: various bone fragment fixation, ligament fixation, non-unions and mal-unions. Please refer to the operative technique for a complete list of indications and contra-indications.

Features & Benefits

- Stardrive recess for 3.0/4.0mm cannulated screw, hexagonal recess for 4.5/7.3mm cannulated screw
- Cannulated screws may be inserted percutaneously over a guide wire (threaded and non-threaded guide wires are available)
- Cancellous thread profile for maximum pull-out resistance, also provides greater inter-fragmentary compression
- Reverse-cutting flutes facilities implant removal
- Each specification has short and long threaded length (except 4.5mm cannulated screw)
- Self-drilling and self-tapping flutes facilitates screw insertion
- Washer prevents screw head from sinking into osteopenic bone.



Implants

3.0 Cannulated Screw(Threaded Length = 1/3 Screw Length)

P/N(TA)	Screw Length	Threaded Length
44309016E	16mm	5mm
44309017E	17mm	6mm
44309018E	18mm	6mm
44309019E	19mm	6mm
44309020E	20mm	7mm
44309021E	21mm	7mm
44309022E	22mm	7mm
44309023E	23mm	8mm
44309024E	24mm	8mm
44309025E	25mm	8mm
44309026E	26mm	9mm
44309027E	27mm	9mm
44309028E	28mm	9mm
44309029E	29mm	10mm
44309030E	30mm	10mm
44309032E	32mm	11mm
44309034E	34mm	11mm
44309036E	36mm	12mm
44309038E	38mm	12mm
44309040E	40mm	12mm

3.0 Cannulated Screw(Threaded Length = 1/2 Screw Length)

P/N(TA)	Screw Length	Threaded Length
44310016E	16mm	8mm
44310017E	17mm	8mm
44310018E	18mm	9mm
44310019E	19mm	9mm
44310020E	20mm	10mm
44310021E	21mm	10mm
44310022E	22mm	11mm
44310023E	23mm	11mm
44310024E	24mm	12mm
44310025E	25mm	12mm
44310026E	26mm	13mm
44310027E	27mm	13mm
44310028E	28mm	14mm
44310029E	29mm	14mm
44310030E	30mm	15mm
44310032E	32mm	15mm
44310034E	34mm	17mm
44310036E	36mm	17mm
44310038E	38mm	19mm
44310040E	40mm	19mm

4.0 Cannulated Screw(Threaded Length = 1/3 Screw Length)

P/N(TA)	Screw Length	Threaded Length
44307022E	22mm	8mm
44307024E	24mm	8mm
44307026E	26mm	9mm
44307028E	28mm	9mm
44307030E	30mm	10mm
44307032E	32mm	10mm
44307034E	34mm	12mm
44307036E	36mm	12mm
44307038E	38mm	13mm
44307040E	40mm	13mm
44307042E	42mm	14mm
44307044E	44mm	14mm
44307046E	46mm	16mm
44307048E	48mm	16mm
44307050E	50mm	17mm
44307052E	52mm	17mm

4.0 Cannulated Screw(Threaded Length = 1/3 Screw Length)

P/N(TA)	Screw Length	Threaded Length
44308022E	22mm	11mm
44308024E	24mm	12mm
44308026E	26mm	13mm
44308028E	28mm	14mm
44308030E	30mm	15mm
44308032E	32mm	16mm
44308034E	34mm	17mm
44308036E	36mm	18mm
44308038E	38mm	19mm
44308040E	40mm	20mm
44308042E	42mm	21mm
44308044E	44mm	22mm
44308046E	46mm	23mm
44308048E	48mm	24mm
44308050E	50mm	25mm
44308052E	52mm	26mm

4.5 Cannulated Screw

P/N(TA)	Screw Length
44301022E	22mm
44301024E	24mm
44301026E	26mm
44301028E	28mm
44301030E	30mm
44301032E	32mm
44301034E	34mm
44301036E	36mm
44301038E	38mm
44301040E	40mm
44301042E	42mm
44301044E	44mm
44301046E	46mm
44301048E	48mm
44301050E	50mm
44301052E	52mm
44301054E	54mm
44301056E	56mm
44301058E	58mm
44301060E	60mm

7.3 Cannulated Screw
(Threaded Length 16mm)

P/N(TA)	Screw Length	Threaded Length
44304050E	50mm	16mm
44304055E	55mm	16mm
44304060E	60mm	16mm
44304065E	65mm	16mm
44304070E	70mm	16mm
44304075E	75mm	16mm
44304080E	80mm	16mm
44304085E	85mm	16mm
44304090E	90mm	16mm
44304095E	95mm	16mm
44304100E	100mm	16mm
44304105E	105mm	16mm
44304110E	110mm	16mm
44304115E	115mm	16mm

7.3 Cannulated Screw
(Threaded Length 32mm)

P/N(TA)	Screw Length	Threaded Length
44306060E	60mm	32mm
44306065E	65mm	32mm
44306070E	70mm	32mm
44306075E	75mm	32mm
44306080E	80mm	32mm
44306085E	85mm	32mm
44306090E	90mm	32mm
44306095E	95mm	32mm
44306100E	100mm	32mm
44306105E	105mm	32mm
44306110E	110mm	32mm

Washer

P/N(TA)	Specification
44305065E	Ø6.5
44305070E	Ø7
44305010E	Ø10
44305013E	Ø13